7.2 NETTOYAGE

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7.2 CLEANING

Note: before cleaning the bath, disconnect it from the main supply and let the heating element cool down.

- ◆ Empty the tank using the emptying pump (except model 601/3, 602/3 and 602/5).
- ◆ For cleaning the stainless steel tank use a cloth, never use scourers since they could scratch the tank. Take special care when cleaning the heating element and the temperature sensor.
- ◆ To avoid the appearance of calcareous deposits, it is recommended to clean the tank periodically with a mixture of water and vinegar (acetic acid) or with a solution of 10% of 28% chlorhydric acid. Rinse with clean water afterwards.
- ◆ Do not use the bath if the heating element is covered with lime since it can lose heating power and deteriorate more quickly. No utilizar el baño si la resistencia está cubierta de cal ya que perderá poder de calefacción y se deteriorará más rápidamente.
- ◆ For cleaning the rest of the equipment, it is recommended to use soapy water without abrasive products.

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

8. STORAGE AND TRANSPORT

- ◆ Keep the equipment in a clean and dry environment with good ventilation and free from corrosive gases.
- ◆ During transport, prevent the equipment from wetting and crashing.

**INSTRUCTIONS ABOUT ENVIRONMENT PROTECTION**

Do not throw this equipment in the garbage when its life cycle had been completed; take it to a collection point for recycling electrical and electronic equipments. It does not contain dangerous or toxic elements for human health but it may cause damage to the environment.

Materials may be recycled as indicated. If you recycle these materials, you are contributing to environment protection. Please contact the administration of your community to be advised on collection points.

6.4 EMPTYING

All baths, except model 601/3, 602/3 and 602/5, are equipped with a pump and drainpipe (1.4) for an easy emptying.

- ◆ Before emptying the bath make sure the On/Off switch (Power) (1.2) is in off position (O); otherwise the heating element could be working when the tank is empty.
- ◆ Press On/Off switch (Outlet) (1.3) and the pump will start emptying the tank through the drainpipe (1.4).

6.5 AUTO-ADJUST

The auto-adjust function allows the adjustment of the internal parameters of the regulator to a certain working temperature. For this proceed as follows:

- ◆ Select the working temperature as explained on Point 6.1 of this manual.
- ◆ Keep pressed ▼ button (2.3) until the AT indicator (2.9) starts flicking; auto-adjust has begun.
- ◆ Auto-adjust will finish when AT indicator (2.9) lights off.
- ◆ After the auto-adjust the equipment will operate perfectly.

7. MAINTENANCE AND CLEANING

7.1 MAINTENANCE

Proper maintenance helps to keep instruments working in a good state and lengthen its lifetime.

Note: all the processes of use mention above will not have any value unless it is kept a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The equipment is supplied with a standard Schuko wire; it has to be plugged to a an earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.
- ◆ If the equipment is not going to be used during a long period of time; make sure it is unplugged.
- ◆ Never try to repair the equipment by yourself, since warranty will be lost and damages in the equipment, users (burns, hurts,...), electrical installation or nearby equipments can be caused.
- ◆ In the event of breakdown please contact your distributor to overhaul through Nahita Technical Assistance Department.

INDEX OF LANGAGES

| | |
|------------------|---------|
| Castellano | 2 - 10 |
| Inglés | 11 - 19 |
| Francés | 20 - 28 |

INDEX OF CONTENTS

| | |
|-----------------------------------|----|
| Preface and Warranty | 12 |
| 1. Safety insrtuctions | 12 |
| 2. Specifications and use | 13 |
| 3. Inspection | 14 |
| 3.1 Reception | 14 |
| 3.2 Listing of items | 14 |
| 4. Installation | 15 |
| 5. Control panel | 15 |
| 6. Operating mode | 16 |
| 6.1 Temperature setting | 16 |
| 6.2 Time setting | 17 |
| 6.3 Alarm | 17 |
| 6.4 Emptying | 18 |
| 6.5 Auto-adjust | 18 |
| 7. Maintenance and cleaning | 18 |
| 8. Storage and transport | 19 |

PREFACE AND WARRANTY

Thank you for purchasing a Nahita thermostatic water bath. Users should read this manual carefully, follow the instructions and procedures, and beware of all the cautions when using this instrument.

If it is necessary, contact your dealer to get the necessary help and advises for a proper use of the equipments. For this, user must give information about the serial number of the equipment, the problems found, the procedures followed to try to solve the problem and contact information.

Nahita thermostatic water baths are warranted to be free from defects in materials and workmanship under normal use and service, for a period of 24 months from the date of invoice. It shall not apply to any product or parts which have been damaged on account of improper installation, improper connections, misuse, accident or abnormal conditions of operation.

For claim under the warranty please contact your distributor in order to send the equipment to Nahita Technical Service Department.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

1. SAFETY INSTRUCTIONS

- ◆ Read the user's manual before operating with the equipment. This equipment must be used only by qualified staff.
- ◆ Connect the device to an earthed power supply to ensure safety of machine and experiment; the voltage indicated on the equipment must correspond to that of the power supply used.
- ◆ The instrument can only be disconnected from the main power supply by pulling out the connector plug, not the cable.
- ◆ This equipment is forbidden to use with inflammable, explosive, poisonous and strong corrosive solutions.
- ◆ Make sure the equipment is placed and installed in a plain, horizontal and stable surface.
- ◆ Respect the minimum safety distances between devices, between the device and the wall and above the assembly (min.200 mm).
- ◆ Keep the equipment away from high magnetic field.
- ◆ The bath must be only used with the tank full of water; if the heating element operates without being covered by water it will be damaged.
- ◆ Be careful when touching the equipment. The metallic tank can reach temperatures up to 100°C (200°C for oil baths). Once the equipment has been turned off it can remain hot for several minutes

- ◆ Use ▲ (2.2) and ▼ (2.3) buttons to select the desired temperature value (show non SV display). Use ◀ (2.4) button to move from one digit to the next.
- ◆ Press Set (2.1). The heating element will start working to reach the set temperature; PV display will show the real temperature of the water, while SV display will show the set temperature.
- ◆ While the heating element is working, RUN indicator (2.7) will be on (green) and OUT indicator (2.8) will flick (Green); once the set temperature is reached then OUT indicator (2.8) will light off.
- ◆ To get an homogeneous temperature in all the volume of water, stir regularly the water in the tank.
- ◆ To assure a constant and homogeneous temperature and prevent water evaporation from the bath, place the lid and the appropriate reducing rings to the size of the recipient to be warmed.

6.2 TIME SETTING

By default, when turning the equipment on the time will be off, so as the equipment will be operating until it is turned off by pressing On/Off switch (1.2). In case it is necessary to work during periods of time, proceed as follows:

- ◆ Press ▲ button (2.1). Command "St" will be displayed on PV display (2.5) and by default "Off" will be displayed on SV display (2.6).
- ◆ Use ▲ (2.2), ▼ (2.3) and ◀ (2.4) buttons to set the working time value; time value can be set up to 9999 min.
- ◆ Press Set (2.1) to save the time value; the heating element will star working.
- ◆ During operation RUN (2.7) and ST (2.10) indicators will be on (green) and OUT indicator (2.8) will flick (green).
- ◆ Time count will start as soon as set temperature is reached.

Note: if working time is long, keep an eye on water level to avoid the heating element working without being covered by water.

- ◆ When working time is finished the equipment and heating element will stop heating (all indicators will be off). No alarm will sound when time finishes.
- ◆ To restart the bath, keep pressed ◀ button (2.4) until RUN indicator (2.7) lights on. A new working temperature and time can be set (time will be again Off by default).


6.3 ALARM

In case real temperature exceeds in more than 10°C the set temperature, the alarm will sound; at the same time real temperature will be displayed on PV display (2.5), command dHAL will be displayed on SV display (2.6) and ALM indicator (2.11) will flick (red).

- ◆ Turn the equipment off by pressing On/Off switch (Power) (1.2) to stop the alarm sound.

Front panel

- 2.1 Set: parameter setting button
- 2.2 ▲: parameter increase button
- 2.3 ▼: parameter decrease button
- 2.4 ◀: digit selection button
- 2.5 PV display
- 2.6 SV display
- 2.7 RUN indicator
- 2.8 OUT indicator: heating element working
- 2.9 AT indicator: autoadjustment
- 2.10 ST indicator: set time
- 2.11 ALM indicator: alarm

6. OPERATING MODE

Note: before operating with the bath, check the level of water in the tank; if it is necessary add more water. Bear in mind that if the heating element is working without being covered by water it will be damaged.

- ◆ Press the On/Off switch (1.2).
- ◆ PV display (2.5) will show real temperature of the bath, while SV display (2.6) will show the set temperature value.
- ◆ After its use do not forget turning off the bath by pressing the On/Off switch (1.2) and unplug it from the power supply.
- ◆ If the equipment is not going to be used during a long period of time, empty the tank with the pump and drainpipe (1.4).

6.1 TEMPERATURE SETTING

- ◆ The maximum working temperature is 100° C (200 °C for oil baths); proceed as follows to set working temperature:
- ◆ Press Set (2.1); command "Sv" will be displayed on PV display (2.5) and temperature set value on SV display (2.6).

- ◆ Pay attention to the set temperature; never leave unattended the equipment during operation and with heating function on.
- ◆ Temperature must always be set to at least 50°C lower than the ignition point of the media used.
- ◆ Pay attention to the safety and resistance to high temperatures of the container used. Process pathogenic materials only in closed vessels.
- ◆ During operation wear the necessary personal protections to avoid the risk of possible damages as:
 - Burns caused by splashing and evaporation of liquids
 - Intoxication caused by release of toxic or combustible gases.
- ◆ Before operation, always check the instrument and accessories are not damaged. Do not use damaged components. Safe operation is only guaranteed with the accessories described in the "Listing of items" chapter. Accessories must be securely attached to the device and cannot come off by themselves. Always disconnect the plug before fitting accessories.

2. SPECIFICATIONS AND USE

- ◆ Nahita digital thermostatic water baths have been designed to get a homogeneous temperature in samples and solutions.
- ◆ For a better performance, Nahita baths are thermally isolated between the tank and the external case, allowing a higher general performance, a better temperature homogenization and energy saving. The depth of the tank allows a bigger immersion and contact of samples with heating source, thus being appropriated to use in Biology, Chemistry, etc.
- ◆ Thermostatic baths present an oven-painted external case that contains a stainless steel tank made out of one single piece of welds; inside this tank, at the lower part, there are the temperature sensor and the heating element, covered by a metallic rack to avoid direct contact between the recipients and the heating element.
- ◆ The tank can be covered by the lid with holes and thanks to the reducing rings, get the appropriate diameter of the hole for any type of recipient. The diameter of the holes is approx. 11 cm and diameter of reducing rings is 8 and 5 cm.
- ◆ Baths are equipped with a pump and drainpipe to be easily emptied (except model 601/3, 602/3 and 602/5).

Technical features chart

| Series 601 | 601/3 | 601/5 | 601/12 | 601/19 |
|--------------------------|----------------------------------|--------------|----------|----------|
| Water baths | | | | |
| Reference | 50601003 | 50601005 | 50601012 | 50601019 |
| Series 602 | 602/3 | 602/5 | | |
| Oil baths | | | | |
| Reference | 50602003 | 50602005 | | |
| Capacity | 2.5 L | 5 L | 12 L | 19 L |
| Useful dimensions (cm) | 15x13x10 | 30x15x10 | 32x30x10 | 50x30x10 |
| External dimensions (cm) | 18x17x21 | 33x18x21 | 36x33x22 | 54x34x22 |
| Holes | 1 | 2 | 4 | 6 |
| Max. temperature | RT+5°C-100°C / 200°C (Oil baths) | | | |
| Resolution | 1 °C | | | |
| Accuracy | ± 1 °C | | | |
| Timer | 1' - 9999' | | | |
| Consumption | 400 W | 500 W | 1000 W | 1500 W |
| Fuse | 3 A | 5 A | 8 A | 10 A |
| Power | 220 V ± 10%; 50/60 ± 1 Hz | | | |

3. INSPECTION

3.1 RECEPTION

- ◆ Unpack the equipment carefully and check for any damages which may have arisen during transport. If it happens, please contact your distributor.
- ◆ If there is any apparent damage on the equipment, do not plug it into the power supply.
- ◆ Keep the original wrapping; you will always need it for returns enclosed with all original accessories. We will only accept the equipment return provided it comes in its original wrapping.

3.2 LISTING OF ITEMS

Digital thermostatic water baths are supplied with the following accessories:

- | | | | |
|------------------|------|---------------------------------|--------------------------|
| - Bath | 1 pc | - Reducing rings with lid | Depending on model |
| - Power cable | 1 pc | - Rack to cover heating element | .. 1pc |
| - Lid with holes | 1 pc | - User's manual | 1 pc |

4. INSTALLATION

- ◆ Place the bath on a horizontal, plane and stable surface appropriate for the weight of the equipment.
- ◆ Make sure the switches of the equipment are in Off (0) position.
- ◆ Make sure the required operating voltage and power supply voltage match; ensure the socket is earthed reliably.
- ◆ Put the metallic rack so as it covers the heating element and fill the tank up to 2/3 its capacity.
- ◆ It is recommended to use decalcified water.

Note: DO NOT use deionized or demineralized water since they have a high rusting power.

- ◆ To prolong bath lifetime and avoid rusting, dissolve 1 g of CO₃Na₂ (sodium carbonate) per liter of water used in the bath.

5. CONTROL PANEL

- 1.1 Control panel
- 1.2 On/Off (Power): general on/off switch
- 1.3 On/Off (Outlet): emptying pump on/off switch
- 1.4 Drainpipe
- 1.5 Socket