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Hersteller

Symbol description

WARNING: Failure to observe may result in injury or even death.

CAUTION: Failure to observe may result in damage to or even destruction of the product.

NOTE: Special information on the operation of the instrument.

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1 Intended use

The SmartMex is intended to be used as a suture retention device to facilitate percutaneous or endoscopic soft tissue procedures.

1.1 Indications

- Meniscal repair and allograft transplantation procedures
- Anchoring the allograft to the meniscal rim during meniscal allograft transplantation procedures.

1.2 Contraindications

- Pathological conditions in the soft tissue that would prevent secure fixation of the device.
- Known hypersensitivity to the implant material. Where material sensitivity is suspected, appropriate tests should be carried out and sensitivity ruled out prior to implantation.
- Do not use for surgeries other than those indicated.

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2 Factors which may impair the success of the operation

- Inadequate soft tissue quality (e.g. previous operation)
- Systemic diseases or metabolic disorders
- History of infections
- Drug dependence and abuse of alcohol and medicaments
- Mental incapacity of patient to understand the instructions of the physician and to comply with them
- Implant migration

Possible adverse effects

- I. Infection, both deep and superficial
- II. Allergies, mild inflammatory and foreign body reactions to implant material
- III. Anchorage slippage
- IV. Anchor loosening
- V. Broken sutures
- VI. Lesion of soft tissue/nerve
- VII. Injury to cartilage
- VIII. Pain

3 Patient information

The physician must inform the patient of the factors which may compromise the success of the operation and possible adverse effects as listed under items “Contraindications” and “Possible adverse effects” as far as they apply to the individual patient. The patient should also be informed about any measures he can take to minimize said factors and adverse effects and the possible consequences.

4 Surgical technique

Surgeons who are not familiar with the specific product and instrument set are expected to study the surgical technique thoroughly before performing the operation. The operation technique is available on the KARL STORZ Homepage under the following number:

- 96156620DF

5 Postoperative strategy

The device is only used to assist the normal healing process. Special weight-bearing and flexion limits for the joint may be applied. The patient must be instructed accordingly.

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6 MRI safety information

Non-clinical testing of the product has demonstrated that the implant is MR safe.

A patient with this implant can be safely scanned in an MR system under any conditions.

In non-clinical testing, the image artefacts caused by the tested devices extend 1 mm or less for a gradient echo and spin echo sequence as provided in the test standard.

7 Safety information

 **NOTE:** This product is packaged in protective packaging and is sterile.

 **WARNING:** Risk of injury and damage to the products: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals and always follow the instructions given precisely. Check the functioning of the products used in combination.

 **WARNING:** Risk of injury: Incorrect application of medical instruments poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

 **WARNING:** Risk of injury: The product has to be used only by experienced surgeons with specific training in the use of this device because this is a technically demanding procedure presenting a risk of serious injury to the patient.

 **WARNING:** Before use, check the expiration date and ensure that the packaging is undamaged. Products which have expired or which are damaged must not be used.

 **WARNING:** Risk of infection: The use of non-sterile implants poses a risk of infection for patients. Do not remove the products from the sterile packaging until immediately before use.

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WARNING: Risk of infection: Reprocessing of disposable articles can jeopardize the functionality and safety of the product and lead to infections. Disposable articles must not be reprocessed. Dispose of disposable articles in accordance with the applicable local regulations.

WARNING: Risk of injury: Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional raw surfaces, sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.

WARNING: Risk of injury: Overloading the instrument by exerting too much force on it may cause the medical device to break, bend, and malfunction, and consequently injure the patient or user. Do not overload the instruments. Do not bend bent instruments back to their original positions.

NOTE: Any decision to remove the implant should take into account potential risk associated with a second surgical procedure. Implant removal should be followed by adequate postoperative management.

NOTE: Until healing is complete, the fixation should be considered as temporary and may not support weight-bearing or other unsupported stresses.

NOTE: The patient should be fully informed about the need to limit activities during healing.

NOTE: Care should be taken to avoid damage when handling the suture. Avoid crushing or crimping the suture strand that might occur as the result of grasping the fasteners with instruments.

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8 Application

It is recommended that the suitability of the products for the intended procedure be checked prior to use.

9 Traceability

The self-adhesive identification labels enclosed with the original packaging must be kept so as to ensure traceability and must be attached to the operation report and the patient dossier.

The person responsible for the hospital's material store must keep a record of the devices in stock, including the article no. (REF) and lot no. (LOT) so that rapid action can be taken in the event of changes or recalls. The manufacturer recommends complete clinical, radiological and statistical documentation. The LOT no. and implant type and material must be recorded in the patient file, operation report and implant card (enclosed adhesive labels).

10 Packaging and sterility

Each individual product is single packed and sterilized using gas (EO). Sterility is guaranteed up to the expiry date, provided that the packaging has not been opened or damaged.

11 Product material

The SmartMex is made of the following components:

- Implantable components

Component	Material
Anchor	VESTAKEEP® i4 G
Suture	Dyneema Purity®

- Further disposable components

Component	Material
Sheet	Stainless steel
Needle	Stainless steel
Protection Sleeve	VESTAMID-Care-ME
Handle/Trigger	CYROLITE G-20

NOTE: It is a sterile disposable product.

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13 Storage conditions

The product must be stored in the original unopened packaging in a dry place at room temperature.

14 Liability

Complications or other adverse effects resulting from incorrect indication or surgical procedure, unsuitable choice or handling of products or instruments, lack of sterility at the site of operation, etc. are the responsibility of the surgeon.

15 Disposal

National laws and regulations and the policies of the clinic to dispose of medical waste must be observed.

16 Product specification and reference numbers

The article number of the SmartMex, All-Inside Meniscal Repair System is:

- 28729M – SmartMex, All-Inside Meniscal Repair System

The following instruments are recommended for use with the product:

- 28728SH – Half Pipe
- 2870356KP – Knot Pusher and Suture Cutter

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CE 0123

17 Handling the accessories

WARNING: Risk of infection: These instruments are not sterile upon delivery. The use of non-sterile instruments poses a risk of infection for patients, users and third parties. Inspect instruments for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the instruments before initial use and before and after every subsequent use using validated procedures.

WARNING: Risk of infection: Incorrectly reprocessed instruments expose patients, users and third parties to a risk of infection as well as the danger of an instrument malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

NOTE: The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' (Art. no. 96216003) can be downloaded or requested by visiting www.karlstorz.com.

WARNING: Risk of injury: Overloading the instrument by exerting too much force may cause the medical device to break, bend, and malfunction, and consequently injure the patient or user. Do not overload the instruments. Do not bend instruments back to their original positions.

18 Directive compliance

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

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19 Packaging symbols

	Manufacturer
	Date of manufacture
	Use-by date
	Batch code
	Catalogue number
	Number of products in the product packaging
	Sterilized using ethylene oxide
	Do not resterilize
	Do not reuse
	Do not use if package is damaged
	Keep away from sunlight
	Keep dry
	Consult instructions for use
	Federal (USA) law restricts this device to sale by or on the order of a physician.
	MR safe

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