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T Cell Activation Bioassay (TCR $\alpha\beta$ -KO), Propagation Model

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1. United States Food and Drug Administration “Approved Cellular and Gene Therapy Products” web page: www.fda.gov/vaccines-blood-biologics/cellular-gene-therapy-products/approved-cellular-and-gene-therapy-products (Internet; accessed April 2, 2021).
2. van Loenen, M.M. *et al.* (2010) Mixed T cell receptor dimers harbor potentially harmful neoreactivity. *Proc. Natl. Acad. Sci. USA* **107**, 10972-7.
3. Ahmadi, M. *et al.* (2011) CD3 limits the efficacy of TCR gene therapy in vivo. *Blood* **118**, 3528-37.
4. Legut, M. *et al.* (2018) CRISPR-mediated TCR replacement generates superior anticancer transgenic T cells. *Blood* **131**, 311-22.
5. Hebeisen, M. *et al.* (2013) Molecular insights for optimizing T cell receptor specificity against cancer. *Front. Immunol.* **4**, 154.
6. Zhong, S. *et al.* (2013) T-cell receptor affinity and avidity defines antitumor responses and autoimmunity in T-cell immunotherapy. *Proc. Natl. Acad. Sci. USA* **110**, 6973-8.
7. Labrecque, N. *et al.* (2001) How much TCR does a T cell need? *Immunity* **15**, 71-82.
8. Thomas, S. *et al.* (2010) Molecular immunology lessons from therapeutic T-cell receptor gene transfer. *Immunology* **129**, 170-7.

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Fc Effector Bioassays

产品	规格	目录号
ADCC Reporter Bioassay, Complete Kit (Raji)*	1 each	G7015
ADCC Reporter Bioassay, Core Kit*	1 each	G7010
ADCC Reporter Bioassay, F Variant, Core Kit**	1 each	G9790
ADCC Reporter Bioassay, Target Kit (Raji)*	1 each	G7016
FcγRIIIa-H ADCP Reporter Bioassay, Complete Kit**	1 each	G9901
FcγRIIIa-H ADCP Reporter Bioassay, Core Kit**	1 each	G9991
Mouse FcγRIV ADCC Bioassay, Complete Kit	1 each	M1201
Mouse FcγRIV ADCC Bioassay, Core Kit	1 each	M1211

*For Research Use Only. Not for use in diagnostic procedures.

**Not for Medical Diagnostic Use.

Additional kit formats are available.

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Not for Medical Diagnostic Use. Additional kit formats and sizes are available.

Immune Checkpoint Bioassays

产品	规格	目录号
4-1BB Bioassay	1 each	JA2351
CD28 Bioassay	1 each	JA6701
CD28 Blockade Bioassay	1 each	JA6101
CD40 Bioassay	1 each	JA2151
CTLA-4 Blockade Bioassay	1 each	JA3001
GITR Bioassay	1 each	JA2291
ICOS Bioassay	1 each	JA6801
ICOS Blockade Bioassay	1 each	JA6001
LAG-3/MHCII Blockade Bioassay	1 each	JA1111
OX40 Bioassay	1 each	JA2191
PD-1/PD-L1 Blockade Bioassay	1 each	J1250
PD-1+TIGIT Combination Bioassay	1 each	J2211
PD-L1 Negative Cells	1 each	J1191
TIGIT/CD155 Blockade Bioassay	1 each	J2201

Not for Medical Diagnostic Use. Additional kit formats are available.

T Cell Activation Bioassays

产品	规格	目录号
T Cell Activation Bioassay (IL-2)	1 each	J1651
T Cell Activation Bioassay (NFAT)	1 each	J1621

Not for Medical Diagnostic Use. Additional kit formats are available.

Cytokine and Growth Factor Bioassays

产品	规格	目录号
IL-2 Bioassay	1 each	JA2201
IL-6 Bioassay	1 each	JA2501
IL-12 Bioassay	1 each	JA2601
IL-15 Bioassay	1 each	JA2011
IL-23 Bioassay	1 each	JA2511
RANKL Bioassay	1 each	JA2701
VEGF Bioassay	1 each	GA2001

Not for Medical Diagnostic Use. Additional kit formats are available.

Control Antibodies and Proteins

产品	规格	目录号
Control Ab, Anti-4-1BB	50µg	K1161
Control Ab, Anti-CD20	5 µg	GA1130
Control Ab, Anti-OX40	50µg	K1191
Control Ab, Anti-CD40	50µg	K1181
Control Ab, Anti-CTLA-4	100µg	JA1020
Control Ab, Anti-LAG-3	100µg	K1150
Control Ab, Anti-PD-1	100µg	J1201
Control Ab, Anti-TIGIT	100µg	J2051
Control Ab, Anti-TIM-3	100µg	K1210
Recombinant VEGF	10µg	J2371

Detection Reagent

产品	规格	目录号
Bio-Glo™ Luciferase Assay System	10ml	G7941
	100ml	G7940
Bio-Glo-NL™ Luciferase Assay System	10ml	J3081
	100ml	J3082
	1,000ml	J3083

Not for Medical Diagnostic Use.

Luminometers

产品	规格	目录号
GloMax® Navigator System	1 each	GM2000
GloMax® Discover System	1 each	GM3000
GloMax® Explorer System	1 each	GM3500

For Research Use Only. Not for use in diagnostic procedures.

Note: Additional Fc Effector, Immune Checkpoint, T Cell Activation and Cytokine Bioassays are available through Promega Elite Access. To view and order Promega Bioassay products visit: www.promega.com/products/reporter-bioassays / or email: eliteaccess@promega.com

^(a)NOT FOR MEDICAL DIAGNOSTIC USE. FOR IN VITRO USE ONLY. BY USE OF THIS PRODUCT, RECIPIENT AGREES TO BE BOUND BY THE TERMS OF THIS LIMITED USE STATEMENT. If the recipient is not willing to accept the conditions of this limited use statement, and the product is unused, Promega will accept return of the unused product and provide the recipient with a full refund.

This product may not be further sold or transferred by the recipient and may be used only by the recipient, and then only for (1) research use; (2) discovery, development and monitoring of biologic drugs and vaccines; (3) quality assurance testing of biologic drugs and vaccines; and (4) product release assays for biologic drugs and vaccines. No other commercial use is allowed. "Commercial use" means any and all uses of this product by recipient for monetary or other consideration, including providing a service, information or data to unaffiliated third parties, and resale of this product for any use. Recipient may genetically engineer the cell line using exogenous nucleic acid sequences provided by recipient. Recipient may propagate and store the cells for long-term use. In addition, recipient must use Bio-Glo-NL™ Luciferase Assay System purchased from Promega for all luminescence assays using this product or contact Promega to obtain a license for use of this product with reagents other than Promega's. PROMEGA MAKES NO REPRESENTATIONS OR WARRANTIES OF ANY KIND, EITHER EXPRESSED OR IMPLIED, INCLUDING AS TO MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE WITH REGARD TO THIS PRODUCT. The terms of this agreement shall be governed under the laws of the State of Wisconsin, USA.

^(b)Broad Limited License.

By purchase of this material, Promega Corporation ("Promega") conveys to purchaser (the "Limited Licensee") the non-transferable right to use the material purchased from Promega and derivatives made by Limited Licensee solely for research conducted by Limited Licensee in accordance with the following requirements: (i) Limited Licensee shall not sell or otherwise transfer material or derivatives to any other person or entity, or use material or derivatives to perform services for the benefit of any other person or entity, except to transfer the material or derivatives to a Limited Licensee affiliate provided such affiliate agrees to be bound by the terms of this Limited License; (ii) Limited Licensee shall use material and derivatives only for its internal research within the Customer Field, which may include internal research within the Customer Field in connection with product research, and not for any Commercial Purposes; (iii) Limited Licensee shall use material and derivatives in compliance with all applicable laws and regulations, including without limitation applicable human health and animal welfare laws and regulations; (iv) Institutions shall provide no warranties of any kind to Limited Licensee (statutory or implied) concerning the Intellectual Property or material or derivatives, including without limitation, as to product quality, condition, description, merchantability, fitness for a particular purpose, noninfringement of intellectual property rights or the absence of latent or other defects, and all such warranties are hereby expressly disclaimed; (v) Institutions shall expressly disclaim any warranty regarding results obtained through the use of material or derivatives, including without limitation any claim of inaccurate, invalid or incomplete results; (vi) Institutions and their directors, trustees, officers, employees, agents, faculty, affiliated investigators and students shall have no liability to Limited Licensee, including without limitation, for any loss of use or profits, business interruption or any consequential, incidental, special or other indirect damages of any kind, regardless of how caused and regardless of whether an action in contract, tort, strict product liability or otherwise; (vii) Limited Licensee shall indemnify, defend and hold harmless the Indemnitees and HHMI Indemnitees against any liability, damage, loss or expense (including without limitation reasonable attorneys' fees and expenses) incurred by or imposed upon any of the Indemnitees or HHMI Indemnitees, as applicable, in connection with any claims, suits, investigations, actions, demands or judgments arising out of or related to the exercise of any rights granted to Limited Licensee under this limited license or any breach of limited license by such Limited Licensee, provided that, to the extent the foregoing is not permitted by law, Limited Licensee agrees, to the extent permitted by law, that it, and not the Indemnitees or HHMI Indemnitees, as applicable, shall be responsible for any liability, damage, loss or expense arising out of or related to the exercise of any rights granted to Limited Licensee under this limited license or any breach of this limited license by Limited Licensee; and (viii) material and derivatives and their use may be the subject of one or more issued patents or pending patent applications owned by one or more Institutions, and the purchase of material does not convey a license under any claims in the foregoing patents or patent applications directed to the material or derivatives or use, production or commercialization thereof, except as expressly set forth in this limited license. Nothing in this limited license shall be construed to confer any rights upon Limited Licensee by implication, estoppel or otherwise as to any technology or patent rights of Broad or any other entity other than the Intellectual Property. In addition, nothing in this limited license shall be construed to confer upon Limited Licensee or any third party any rights under or to the Intellectual Property outside of the Customer Field.

Commercial Purposes means (a) the practice, performance or provision of any methods, process or service; or (b) the manufacture, production, sale, use, distribution, disposition or importing of any product, in each case (a) or (b) for consideration of any kind, for the purpose of sale or commercial exploitation, or on any other commercial basis.

Customer Field means use as a research tool for research purposes and use for quality control/quality assurance testing and product lot release assays, provided, however, that notwithstanding the foregoing, the Customer Field shall expressly exclude:

- 1) Any human or clinical use, including, without limitation, any administration into humans or any diagnostic or prognostic use;
- 2) Any human germline modification, including modifying the DNA of human embryos or human reproductive cells;
- 3) Any in vivo veterinary or livestock use;
- 4) The development, manufacture (including any bioproduction), distribution, importation, exportation, transportation, sale, offer for sale, marketing, promotion or other exploitation or use of the Intellectual Property or the material or derivatives as a therapeutic or diagnostic for humans or animals;
- 5) Products that provide nutritional benefits and are regulated by a regulatory authority as a drug or biologic pursuant to Section 505 of the Federal Food, Drug, and Cosmetic Act of 1938, as amended, Section 351 of the Public Health Service Act of 1944, as amended, or any successor laws, or equivalent laws or regulations in jurisdiction outside the United States;
- 6) Any agricultural use, including but not limited to, the use or application in the cultivation, growth, manufacture, exportation, or production of any tobacco product; and
- 7) Any use or application relating to gene drive.

HHMI Indemnitees means Howard Hughes Medical Institute and its trustees, officers, employees and agents.

Indemnitees means each Institution, their affiliates, and their current and former trustees, directors, officers, faculty, affiliated investigators, students, employees, medical and professional staff and agents and their respective successors, heirs and assigns.

Institutions means each of The Broad Institute, Inc., Massachusetts Institute of Technology, Harvard College, University of Iowa Research Foundation, University of Tokyo, The Rockefeller University, New York University, New York Genome Center, and Whitehead Institute of Biomedical Research collectively.

Intellectual Property means (a) any patents and patent applications listed on Appendix A; (b) any substitutions, divisions, continuations, reissues, renewals, re-examinations or extensions of any patents or patent applications of (a), any continuations-in-part of any patents or patent applications of (a) to the extent the claims are directed to subject matter specifically described in a patent or patent application set forth on Appendix A, and any application claiming priority to a provisional application set forth on Appendix A to the extent the claims are directed to subject matter specifically described in a provisional patent application set forth on Appendix A; (c) any letters patents and/or the equivalent issuing therefrom in any jurisdiction based on, and only to the extent of, any of the foregoing; and (d) any foreign or international equivalents of any of the foregoing, including without limitation any applications filed under the Patent Cooperation Treaty.

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* Broad does not purport to grant any rights in the asterisked applicant's interest in these applications in this Agreement.

^(c) This Product and its use are the subject of one or more of the following patents and applications controlled by Sigma-Aldrich Co., LLC (SIGMA): Patent applications and issued patents that entered their respective National Stage from PCT International Pub. No. WO 2014/089290, including, but not limited to, the following, and substitutions, divisions, continuations, continuations-in-part, reissues, reexaminations, and extensions thereof: Australia Patent Nos. 2013355214; 2017204031; and 2018229489; Canada Patent No. 2,891,347; China Patent No. CN105142669; European Patent Nos. EP 3 138 910 B1, 3 138 911 B1, and EP 3 138 910 B1; Israel Patent No. IL238856; Singapore Patent No. 11201503824S; South Korea Patent No. 10-1844123; and U.S. Patent Application Serial Nos. 15/188,911; 15/188,924; 15/188,927; 15/188,931; and 15/456,204 (the "Patent Rights").

BEFORE OPENING OR USING THIS PRODUCT, PLEASE READ THE TERMS AND CONDITIONS SET FORTH IN THIS LICENSE AGREEMENT. YOUR USE OF THIS PRODUCT SHALL CONSTITUTE ACKNOWLEDGMENT AND ACCEPTANCE OF THESE TERMS AND CONDITIONS. If you do not agree to use this Product pursuant to the terms and conditions set out in this License Agreement, please contact Promega Corporation within ten (10) days of receipt to return the unused and unopened Product for a full refund, provided, however, that custom-made Products may not be returned for a refund.

The purchase of this Product conveys to you (the "Buyer") the NON-TRANSFERABLE right to use the Product for Licensed Research Use (see definition below) subject to the conditions set out in this License Agreement. If you wish to use this Product for any purpose other than Licensed Research Use, you must first obtain an appropriate license (see information set out below).

This Product may not be used for any purpose other than Licensed Research Use. Your right to use this Product for Licensed Research Use is subject to the following conditions and restrictions:

1. "Licensed Research Use" means any use for research purposes, with the following exceptions: (i) Buyer may not sell or otherwise transfer the Product (including without limitation any material that contains the Product in whole or part) or any Related Material to any other third party; (ii) Buyer may use the Product, components of the Product, and any Related Material, only for your internal research within the Field, and not for any Commercial Purposes (for any use of Product, components of a Product and any Related Material by any entity, including but not limited to contract service providers ("CRO") offering services for consideration, please contact Promega for additional information), and this License Agreement does not grant or in any way convey any rights to Buyer to further edit Products, components of Products or Related Materials under any claims in the foregoing patents and patent applications; (iii) Buyer shall use the Product and any Related Material in compliance with all applicable laws and regulations, including without limitation applicable human health and animal welfare laws and regulations; (iv) SIGMA shall provide no warranties of any kind to the Buyer (statutory or implied) concerning the Patent Rights, the Product or any Related Materials, including without limitation, as to product quality, condition, description, merchantability, fitness for a particular purpose, noninfringement of intellectual property rights or the absence of latent or other defects, and all such warranties are hereby expressly disclaimed; (v) SIGMA shall expressly disclaim any warranty regarding results obtained through the use of the Product or Related Materials, including without limitation any claim of inaccurate, invalid or incomplete results; (vi) SIGMA, and their directors, officers, employees and agents, shall have no liability to the Buyer, including, without limitation, for any loss of use or profits, business interruption or any consequential, incidental, special or other indirect damages of any kind, regardless of how caused and regardless of whether an action in contract, tort, strict product liability or otherwise; (vii) the Buyer shall indemnify, defend, and hold harmless SIGMA and their current and former directors, officers, employees and agents, and their respective successors, heirs and assigns (the "Indemnitees") against any liability, damage, loss, or expense (including without limitation reasonable attorneys' fees and expenses) incurred by or imposed upon any of the Indemnitees in connection with any claims, suits, investigations, actions, demands or judgments arising out of or related to the exercise of any rights granted to the Buyer hereunder or any breach of this License Agreement by such Buyer, provided that, to the extent the foregoing is not permitted by law, the Buyer agrees, to the extent permitted by law, that it, and not the Indemnitees shall be responsible for any liability, damage, loss or expense arising out of or related to the exercise of any rights granted to the Buyer under this License Agreement or any breach of this License Agreement by Buyer, and (viii) the Product and its creation and/or use may be the subject of one or more issued patents and/or pending patent applications owned by SIGMA and the purchase of the Product and its use as defined in this License Agreement does not convey to the Buyer any license under any claims in the foregoing patents or patent applications except for the limited rights as described in this License Agreement directed toward use of the Licensed Product or commercialization thereof.

2. For purposes of Section 1 above, the following definitions shall apply:

- "Commercial Purposes" means (a) the practice, performance or provision of any method, process or service or (b) the manufacture, sale, distribution, disposition or importing of any product, in each case (a) or (b) for consideration, e.g., a fee, or on any other commercial basis.
- "Field" means use as a research tool for research purposes which includes, but is not limited to drug discovery and development, quality control/ quality assurance testing and product lot release assays, provided, however, that notwithstanding the foregoing, the Field shall expressly exclude (a) any in vivo and ex vivo human or clinical use, including, without limitation, any administration into humans or any diagnostic or prognostic use, (b) the creation of transgenic rodent models and/or derivatives thereof (such derivatives including, but not limited to, rodents' cells and rodents' organs), (c) any in vivo veterinary or livestock use, or non-research agricultural use, or (d) the manufacture, distribution, importation, exportation, transportation, sale, offer for sale, marketing, promotion or other exploitation or use of, or as, a therapeutic or diagnostic for humans or animals.

- “Related Materials” means any progeny, modification or derivative of a Product.

3. Your right to use the Product will terminate immediately if you fail to comply with these terms and conditions. For information on purchasing a license (a) to this Product for purposes other than Licensed Research Use, or (b) to the foregoing patents and patent applications generally, contact your local SIGMA sales representative, who will refer you to the proper licensing representative, or in the USA call 800-325-3010.

^(d)BEFORE OPENING OR USING THIS PRODUCT, PLEASE READ THE TERMS AND CONDITIONS SET FORTH IN THIS LIMITED USE LABEL LICENSE. YOUR USE OF THIS PRODUCT SHALL CONSTITUTE ACKNOWLEDGMENT AND ACCEPTANCE OF THESE TERMS AND CONDITIONS. If you do not agree to use this Product pursuant to the terms and conditions set out in this Limited Use Label License, please contact Promega Corporation within ten (10) calendar days of receipt to return the unused and unopened Product for a full refund; provided, however, that custom-made Products may not be returned for a refund.

Promega is selling this Product under a research use only license from Caribou Biosciences, Inc., under certain patents and patent applications owned or controlled by Caribou. For information on obtaining rights beyond those set forth in this Limited Use Label License please contact licensing@promega.com or Promega Corporation, 2800 Woods Hollow Road, Madison, WI 53711 USA, Attn: Licensing.

The purchase of this Product conveys to you, the Buyer, the NON-TRANSFERABLE right to use this Product solely for Research Use, subject to the conditions set out in this Limited Use Label License. If you wish to use this Product for any purpose other than Research Use, you must first obtain an appropriate license. This Product may not be used for any purpose other than Research Use. Your right to use this Product for Research Use is subject to the following conditions and restrictions:

Research Use means any use for internal research purposes, (e.g., for screening, characterization, quality control/quality assurance testing, product lot release assays, and/or toxicology purposes) for the sole benefit of you, the Buyer.