Specification

Power	AA-size batteries/AC adapter(100-240V/50-60Hz)
Display mode	Digital LCD
Measurement mode	Oscillographic testing mode
Measurement range	Pressure: 0kPa~37.3kPa(0mmHg~280mmHg) Pulse value: (30~160)times/minute
Accuracy	Pressure: 15℃~25℃ within ±0.4kPa(3mmHg) 5℃~40℃ (out of 15℃~25℃) within ±0.8kPa(6mmHg) Pulse value: within ±5%
Speed of pressure increasing	0.5kPa/s~1.1kPa/s(3.75mmHg/s~8.25mmHg/s)
Normal working condition	Temperature: 5℃~40℃ Relative humidity: ≤80%
Storage & transportation condition	Temperature: -20℃~55℃ Relative humidity: 10%~93%
Measurable perimeter of the upper arm	About 20cm~34cm
Weight	Approx. 240g(Excluding the dry cells)
External dimensions	Approx. 520mm×208mm×520mm
Classification of safety II	Internally powered equipment. Type B applied part
Pressure mode	Pressure regulating mode of septum-mode pump
Fast gas-discharging mode	Fast open gas-discharging valve mode
Pressure testing	Semiconductor pressure sensor
Pulse testing	Semiconductor pressure sensor

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
W90303109101
B1208-0
Printed in China

- 感谢您购买本公司产品。
- Thank you for buying our product(s).
- 使用前请您务必仔细阅读本使用说明书,并且要妥善保管。
- Before operating this device,Please read these instructions completely and save this manual for future use.

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EW3109 Upper Arm Blood Pressure Monitor

- This instrument is an automatic electronic sphygmomanometer for measuring upper arm blood pressure. The values of blood pressure and pulse can be easily measured by using oscillography.
- Applicable to the measurement of human body's blood pressure and pulse without traumatism.
- The wrist circumference that can be measured with this wrist blood pressure monitor is about 20 cm - 34 cm.

Our promise:

Any malfunction occurring within of purchase will be dealt with according to the rules in Repair Card attached to the product.

Packing list

- ① Main unit (with pressure cuff) 1
- ② User manual 1
- ③ Repair card 1
- ④ Certificate 1
- ⑤ Record card for blood pressure measurement 1
- ⑥ AA-size batteries 4
- ⑦ AC adapter 1

Included batteries are for device testing. Please to verify the functions of the device during purchase.

Content

Important Safety Instructions	21
Parts Identification	23
Preparation	24
• Cased Battery 24~25
• When using the AC adapter26
How to Use	27
• Positioning the cuff27
• Instructions for correct measurement28
• Methods of measurement29~30
• Storing measurement values31
• Recalling stored Measurement values32
• After use 33
Other	34
• Cautions During Measurements34
• Q&A35~37
• Care And Maintenance	...37
• Troubleshooting38

IMPORTANT SAFETY INSTRUCTIONS

- * Before using this instrument, please read the "IMPORTANT SAFETY INSTRUCTIONS" and be sure to use the instrument properly.
- * The cautions in this section are intended to instruct you in the proper use of this instrument and prevent danger or injury to yourself or other people.

All of the cautions contain important information regarding safety and must be followed.

WARNING: Failure to follow the caution may result in death or serious injury of people.

CAUTION: Failure to follow the caution may result in injury to people, or in damage to property.

WARNING:

Persons with serious obstructions of blood flow in the arm should consult a doctor before using this instrument. Failure to do so may result in deterioration of physical condition.

CAUTION:

- Do not make judgements yourself regarding the measurement results or self-treat yourself.
- If the product does not work or seems strange, stop use immediately and request inspection and service. Failure to do so could result in the risk of fire, electric shock, or burns.
 - There will be large differences between wrist blood pressure values and upper-arm blood pressure values for persons with periphery circulatory obstruction due to diabetes, liver disease, kidney disease, arteriosclerosis, high blood pressure, etc.
 - Take measurements according to your doctor's instructions and have your doctor do the diagnosis.
 - Take medicine according to your doctor's instructions.
- Do not use for purposes other than the measurement of blood pressure. Doing so may result in an accident or malfunction.
- Do not use this instrument on infants or on people who cannot express their feelings. Doing so may result in problems or an accident.
- Do not disassemble this instrument, or attempt to repair or modify it yourself. Doing so may result in an accident or malfunction.

SAVE THESE INSTRUCTIONS


AC ADAPTER NOTES

- When cleaning, unplug the AC adapter from the AC outlet before cleaning. Also, do not plug in or unplug the AC adapter with wet hands. Doing so may result in electric shock or fire.
- Do not spill water on the main body or AC adapter. Doing so may cause electric shock, short-circuiting, or product damage.
- Do not use the AC adapter if it is damaged or if it plugs into the AC outlet loosely. Doing so may result in electric shock or short-circuiting.
- If the product does not operate properly or seems strange, stop use and unplug the AC adapter immediately, and request inspection and service. Failure to do so could result in the risk of an accident, electric shock, or fire.
- When using the AC adapter, make sure it is fully inserted into the AC outlet. Failure to do so could result in electric shock or short-circuiting.
- When the product is not in use, unplug the AC adapter from the AC outlet. Failure to do so could result in deterioration of the insulation due to dust or moisture, causing a leakage electricity fire.
- When unplugging the AC adapter, do not pull on the cord. Always grasp the AC adapter body and unplug it. Failure to do so could result in the risk of electric shock, short-circuiting, or fire.
- Do not insert any object in the AC adapter or allow debris to collect on it. Doing so could result in the risk of electric shock, short-circuiting, or fire.
- Do not scratch, damage, modify, forcibly bend, pull, or twist the AC adapter cord. In addition, do not place heavy objects on top of the cord or pinch it with objects. Do not wrap the cord around the main body and the AC adapter. Doing any of these things may result in the risk of fire or electric shock.
- Do not step on the AC adapter. Doing so may cause an accident or injury.
- The AC adapter is designed for use only with Panasonic Blood Pressure Monitor, model EW3109. Use with other electrical devices may result in death or serious personal injury.
- ALWAYS unplug the AC adapter from the power outlet before performing maintenance.
- Do not allow any material to drape over the AC adapter. Doing so could result in the risk of fire.
- Do not modify or attempt to repair the product or AC adapter by yourself. Doing so could result in the risk of fire, electric shock, or personal injury.

PARTS IDENTIFICATION

Please check the instrument.


Main Body


Optional Accessory: AC adapter

Modelo No. EW31225


INPUT: 100-240V / 50-60Hz
OUTPUT: 6V

Battery of main body


Storage pouch

Specifications of products and accessories are subject to change without notice.

INSTALLING BATTERIES

* Always use alkaline batteries.

1 Turn the main body over, press lightly on the area indicated by (A) on the battery cover, as shown at right and

remove the cover by sliding it in the direction of the arrow.

2 Check the polarities of the batteries and **install** them correctly.


3 **Close** the battery cover until a click heard.

CAUTION:

- Always follow the cautions printed on the batteries.
- Check the polarities of the batteries and install them correctly.
- Promptly remove exhausted batteries.
- If the product will not be used for a long time, remove the batteries.
- Batteries should be used before their expiration date.

Use after the expiration date may result in injury or staining of the surrounding area due to generation of heat, bursting, or leakage.

Replace batteries when:


- The mark shown above appears during use.
- All displays are lit for a moment after the 开始/结束 button is pressed, and then the display goes blank.
- No display appears even if the 开始/结束 button is pressed.

- Approximately 500 measurements (approximately 5 months when measurements are taken 3 times per day) can be taken with a set of fresh Panasonic alkaline batteries (four AA-size LR6 batteries). (Measurement conditions: Room temperature: 23°C ; Pressurized to 170 mmHg; Upper arm circumference: 30 cm)
- When room temperature is low, battery life may become shorter.
- Batteries are included only for testing purposes and with no guarantee of battery life. (In some areas batteries are not included with the unit.)
- A set of new "Panasonic" Zn-Mn batteries will provide about 140 measurements.
- Included Zn-Mn batteries are for instrument testing. Since these batteries are used to ensure the functions of the instrument at the time of purchase, their lifetime may be less than 140 measurements.
- When replacing batteries, always use four fresh alkaline batteries of the same type from the same manufacturer.
- If the instrument will not be used for a long period of time, remove the batteries. Leaving the batteries installed may result in battery leakage and damage to instrument.
- Do not discard batteries with normal garbage. Take to recovery center.

WHEN USING THE AC ADAPTER


Use only the exclusive AC adapter(EW31225)

- 1 Insert tightly the wire in the body of AC adapter.
- 2 Insert the AC adapter plug into the connector on the right side of the main unit.
- 3 Plug the AC adapter fully into the AC outlet.


Be sure to use only the exclusive AC adapter.
Use of other AC adapters could cause an accident.

POSITIONING THE PRESSURE CUFF

- 1 Insert** the plug into the socket.
 - Insert the pressure cuff plug all the way in.

<When wrapping around left arm>


- Position the tube outlet so that it is on the palm side of the arm in line with the little finger of the hand.

Tape fastener

- Please take care of using tape fastener because there's risk of damaging clothes and fabric products.

<When wrapping around right arm>


- 2 Wrap** around upper arm.
 - Wrap around bare arm.

The blood pressure monitor can be used by wrapping the cuff on either upper arm. However, since it has been said that the difference in blood pressure between the left and right arms may be around 10 mmHg, be sure to always use the same arm for measurements.

- 3 Wrap** the sleeve strap. (do not fasten too tightly)
 - Do not wrap clothes into the sleeve strap (otherwise the accuracy of measurement may be affected).
 - Do not wrap too tightly. Loose enough for one finger insertion is appropriate.
 - Fasten the nylon fastener. Make sure the air tube dose not fall out.

FOR ACCURATE MEASUREMENTS

Take measurements while sitting quietly and relaxing.

- 1 Rest your arm elbow** on a table, etc.

- 2** Position your arm so that the pressure cuff is **at the same height** as your heart.

- 3** Position your palm upward, and **relax your hand.**

- ♥ **Take measurements using the upper arm in the same position at the same time every day.**


- Measurement values will vary according to the height of the pressure cuff.
- **It is better to measure your basic blood pressure (the blood pressure when you get up in the morning without pee feeling). If not available, please measure your blood pressure before eating and before intense physical exercise.**
- During measurement, do not touch the main body or the tube.

- ♥ **Relax and sit in a natural posture when taking measurements.**

- Always sit quietly for approximately 4 to 5 minutes before taking measurements.
- Relax your body, arm, and fingers and do not move their muscles.
- Do not talk while taking measurements.
- Please ensure 4 to 5 minutes of intermission between measurements.
- Do not use mobile phone or any other apparatus that may generate electromagnetic wave near the main unit. Otherwise, it may cause misoperation of the blood pressure monitor.

TAKING A READING

Measurement

- Press** the 开始/结束 button.
 - The pressure cuff will automatically be pressurized and a measurement will be taken.
 - The entire display will light for approximately 1 second. If there is air left in the pressure cuff, the fully lit time will be longer.

The heart mark ♥ will blink when pulse is detected during pressurization.

Since measurements are performed during pressurization, there may be changes in operating sounds or the pressurization rate may change during pressurization. This is normal and is not a malfunction

- When measurement has been completed, the blood pressure value and pulse rate will be displayed.


- Pressure will be automatically released from the pressure cuff.
- The blood pressure values and pulse number are displayed at the same time. Blood pressure values are displayed alternately in mmHg and kPa units.
- If the measured value is in the high blood pressure region, the values will blink (for approximately 8 seconds). (For how to read the blood pressure display, see p 30)
- The (M) mark will blink. (For storing measurement values, see p 31)
- *The pulse rate value is the beats per minute calculated based on the pulse rate during measurement.

- To switch off the blood pressure monitor, **press** the 开始/结束 button.
 - If you forget to switch off the blood pressure monitor, it will switch off automatically after approximately 5 minutes.

Reading the blood pressure display

Whether the systolic and diastolic blood pressure values are in the high blood pressure region or normal blood pressure region are judged based on the WHO/ISH (World Health Organization/International Society of Hypertension) blood pressure classifications, and the blood pressure value blinks if required.


If the measured blood pressure value is in the high blood pressure region, the values will blink for approximately 8 seconds.

- Systolic blood pressure: 18.7kPa (140mmHg) or higher (Only immediately after measurement)
- Diastolic blood pressure: 12.0kPa (90mmHg) or higher (1kPa=7.5mmHg)

If measurements could not be taken properly

The remeasure mark "E" will be displayed.

When taking successive measurements, press the 开始/结束 button to switch off the blood pressure monitor once and stay quiet for 4 to 5 minutes, and then take the next measurement.

Advice:

Often with the second measurement, there is less tension and you feel more relaxed, so the blood pressure value is lower. For people who easily get tense, measuring twice each time is very effective. When two measurements are taken, normally the smaller measurement is stored. When taking two measurements, always stay quiet for 4 to 5 minutes after the first measurement before taking the second measurement.

*1 1999 World Health Organization-International Society of Hypertension Guidelines for the Management of Hypertension

STORING MEASUREMENT VALUES

This product can store up to 90 measurement values.


- Stored data will remain in memory even if batteries become exhausted or batteries are removed.

After measurement, **press** the 记忆/调出 button while the (M) mark is blinking.

Storing is completed.

- Up to 90 measurement values can be stored. When the memory is full and additional values are stored, the oldest value in memory will be deleted and the new value will be stored.
- If measurements could not be taken properly (the E mark is shown instead of the blood pressure value), measurement values cannot be stored.
- * However, even if the E mark is shown for the pulse value, the blood pressure value will be stored.

Number of stored data

To erase all measured values stored in memory:

1. Press the 记忆/调出 button to recall stored measurement values.
2. Press the 记忆/调出 button again and hold it pressed until the measurement values disappear (longer than approximately 3 seconds).

* It is not possible to erase single measurement values.

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To switch off the blood pressure monitor,

press the 开始/结束 button.

- If you forget to switch off the blood pressure monitor, it will switch off automatically after approximately 5 minutes.

RECALLING STORED MEASUREMENT VALUES

Measurements can be recalled after measurement is completed or after the power has been switched off.

If no measurement values are stored, no value will be shown.

1 Press the 记忆/调出 button. The average of all measured values in memory will be shown.

- Avg. will be shown at the top right of the display.
- For example, if 90 measured values are stored, the average of 90 measurements will be shown. (If only 1 measurement is stored, the Avg. mark will not be shown.)

2 Press the 记忆/调出 button again. Each time you press the 记忆/调出 button, the previously stored measurement values will be **shown** in order, starting with the newest data.

Press the 记忆/调出 button lightly for only a short time. If you hold the button pressed for approximately 3 seconds or longer, all data in memory will be erased, so please be careful.

3 To cancel measurement value recall, **press** the 开始/结束 button.

- If you forget to switch off the blood pressure monitor, it will switch off automatically after approximately 30 seconds.

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
AFTER USE

Storage


1 Unplug the plug from the main body.

2 Roll the cuff loosely around the main body.

- Please do not wrap the air tube around the main unit.
- Please do not pull or bend the air tube forcibly. Doing so may cause malfunction.


3 Place in the storage pouch.

- Be careful not to bend the tube.

CAUTIONS DURING MEASUREMENTS

In the following situations, accurate measurements cannot be taken.

Within 1 hour of eating or drinking alcohol

Immediately after drinking coffee or tea, or smoking

Immediately after exercising or taking a bath

In a moving vehicle


Wait at least 20 minutes and take a measurement while staying quiet.

In cold places

Take measurements at ambient temperatures of 10°C to 40°C.

When you need to urinate


Wait several minutes after urinating before taking measurements.

Q & A

Q Why are the blood pressure values measured at home different from those measured at the doctor's office?

A. Blood pressure varies 24 hours a day. In addition, it can change greatly depending on the weather, your mental state, exercise, etc. In particular, there is the phenomenon of "White uniform hypertension" at doctor's offices, where the blood pressure measured in a doctor's office is higher because of stress or worry.


In addition, if the blood pressure measured at home is higher than that measured at the doctor's office, the following points should be considered and care should be observed where necessary:

① **Is the pressure cuff wrapped correctly?**

- Make sure that the pressure cuff is not wrapped too loosely or too tightly.
- Is the pressure cuff wrapped around the upper arm?

For how to wrap the pressure cuff correctly, refer to p.27.

② **Do you feel worried or irritated?**

- Before taking measurements, take 2 or 3 deep breaths and relax to stabilize blood pressure. In addition, it is recommended that you sit quietly for about 4 to 5 minutes before taking measurements

Q & A

Q When should I measure my blood pressure?

A. We recommend that you measure your blood pressure 3 times a day.

First time: **When you get up in the morning.** (This is when your blood pressure is most typical of its value during the night and when the movements of your body have had hardly any effect on it.)


Second time: **Before eating lunch.** (This is when your blood pressure is not affected by food and when it is most typical of the highest level of the day.)

Third time: **When you go to bed at night.** (This data value is roughly typical of the average blood pressure value for the day.)

Get to know your blood pressure pattern by measuring your blood pressure in the morning, at midday and in the late evening, and use the data you obtain to manage your blood pressure.

Q Why are measured blood pressure values different each time I take a measurement?

A.

① Blood pressure varies throughout the day, and can also vary depending on how the pressure cuff was wrapped or your posture. Because of this, measurements should be taken at the same time and under the same conditions every day.


② For persons taking blood pressure medication, blood pressure may change greatly due to the effects of the medication.

③ When taking 2 measurements in succession, be sure to stay quiet for 4 to 5 minutes after the first measurement before taking the next measurement.

Q & A /CARE AND MAINTENANCE

Q Why are the measured values different for the left and right wrist?

A.

It is said that differences in blood pressure values between the left and right arms may be around 10 mmHg. The blood pressure can be measured in either the left or right arm, but all measurements should be taken using the same arm.

Care and maintenance

Do not apply excessive force to the instrument or drop the instrument.

- Doing so may cause damage.

If the instrument is stored at temperatures below the freezing point, do not use it immediately. Leave it in a warm place for at least an hour before use.

- If the instrument is not allowed to warm up, it may not pressurize.

Please keep it away from foreign bodies or erosive liquid.

- Doing so may result in damage.

If the instrument becomes dirty, clean it with a soft cloth moistened with warm water or soapy water. (Do not use alcohol, benzene, or paint thinner.)

- Use of such chemicals may result in cracking or discoloration.

For storage, avoid high temperatures, high humidity, and direct sunlight.

- Storage under such conditions may result in damage.

The pressure cuff cannot be washed.

TROUBLESHOOTING

Display	Status before error	Cause and solution
Remeasure E mark is displayed. 	Pressurization to more than 280 mmHg was performed.	Were measurements taken according to the proper procedure and in the correct posture? See pp.27 and 28.
	Pressure decreased suddenly and the remeasure mark E appeared.	
	Measurement was taken to 20 mmHg and then the remeasure mark E appeared.	Was the pressurization cuff wrapped properly? See p.27.
	The heart mark \heartsuit blinked only a few times or not at all. Pressurization doesn't occur even after waiting several dozen seconds.	Tube may not be connected sufficiently. Check the tube connection. See p.27.

Symptom	Possible cause
Systolic or diastolic blood pressure value is high or low.	<ul style="list-style-type: none"> • Position of pressure cuff is low. • Pressure cuff is not wrapped correctly. • Position of pressure cuff is high. • You moved or talked while measurements were being taken. See pp.27 and 28.
Blood pressure value is abnormally high or low.	<ul style="list-style-type: none"> • Posture is different for each measurement. See p.28.
Value is different from that measured by the doctor.	<ul style="list-style-type: none"> • Blood pressure varies in response to minute changes in your mental state, such as your reaction to a doctor's presence. • Sit quietly for 4 to 5 minutes and then take another measurement. See p.28.
Measurement values are different with each measurement.	
No display appears even when the 开始/结束 Start button is pressed.	<ul style="list-style-type: none"> • Batteries are exhausted. • Batteries are installed with the polarities wrong. See pp.24 and 25.
When the 开始/结束 Start button was pressed, the display became lit and then the display disappeared.	
Damage of the main part	<ul style="list-style-type: none"> • Do not repair it yourself. Please contact our designated service agency for electronic sphygmomanometer to examine and repair the device.

***The warranty of this product is only valid in the PRC. (Not available in Hongkong, Macao and Taiwan for the time being.)**