

NOTICE: This document contains references to Varian. Please note that Varian, Inc. is now part of Agilent Technologies. For more information, go to [www.agilent.com/chem](http://www.agilent.com/chem).



## APPLICATION NOTE

Application 1936 - HPLC

### Fatty acids

#### Fatty Acids

Column: Polaris™ 3μ C18-A, 30 x 4.6mm  
Part nr.: A2001030X046  
Eluent: A: H<sub>2</sub>O + 0.1% HCOOH B: MeCN + 0.1% HCOOH A:B = 5:95  
Flow: 1.5 mL/min  
Detection: UV @ 205 nm  
Sample: 1. linolenic acid, 2. linoleic acid, 3. oleic acid  
Injection Volume 2 μL

Three 18-carbon fatty acids are separated in under 52 seconds using an isocratic eluent.

**Oleic** =  $\text{CH}_3(\text{CH}_2)_7\text{CH}=\text{CH}(\text{CH}_2)_7\text{COOH}$

**Linoleic** =  $\text{CH}_3(\text{CH}_2)_4\text{CH}=\text{CHCH}_2\text{CH}=\text{CH}(\text{CH}_2)_7\text{COOH}$

**Linolenic** =  $\text{CH}_3\text{CH}_2\text{CH}=\text{CHCH}_2\text{CH}=\text{CHCH}_2\text{CH}=\text{CH}(\text{CH}_2)_7\text{COOH}$

