

NOTICE: This document contains references to Varian. Please note that Varian, Inc. is now part of Agilent Technologies. For more information, go to www.agilent.com/chem.



Application Note SI-02374

EPA 525.2: Trace Level Determination of Semi-Volatile Organic Compounds in Drinking Water by Liquid-Solid Extraction and Capillary Column GC/MS Using the 240-MS Ion Trap and V:Results™ GC/MS software

Lily Lew, Anaïs Viven and Ed George
Varian, Inc.

Introduction

EPA Method 525.2, "Determination of Organic Compounds in Drinking Water by Liquid-Solid Extraction and Capillary Column Gas Chromatography/Mass Spectrometry," is one of the most commonly used methods for semi-volatile compound analysis in drinking water. This note provides a basic overview of 25 regulated compounds analyzed on the 240-MS GC/MS system for the analysis of water samples by EPA Method 525.2. The 240-MS provides the ultimate in sensitivity, easily exceeding the required detection limits cited in the method. The system can also be configured in both the internal or external ionization mode for added flexibility.




Figure 1. Varian 240-MS Ion Trap Mass Spectrometer with 431-GC.

Instrumentation

- Varian 240-MS Ion Trap Mass Spectrometer
- Varian 431-GC Gas Chromatograph
- Varian 1177 Split/Splitless Injector with Siltek frit insert (Part No. RT210462145)
- 8400 AutoSampler
- V:Results GC/MS software

GC Conditions

Column: FactorFour™ VF-5ms 30 m × 0.25 mm × 0.25 μm
(Part No. CP8944)

Program: 45 °C for 1.5 min, to 240 °C for 2 min at 20 °C/min,
to 300 °C for 4 min at 40 °C/min

MS Conditions

MS Configuration: Internal ionization
Target TIC: 6000 μs
Scan Range: 45-450 m/z
Max Ion Time: 25000 μs
Emission Current: 20 μA
Manifold Temp: 50 °C
Transfer line Temp: 280 °C
Ion Trap Temp: 220 °C

Results and Discussion

A typical total ion chromatogram (TIC) and extracted ion chromatogram (EIC) are shown in Figures 2 and 3.




Figure 2. TIC at 1 ppm of regulated 525.2 mixture, 240-MS internal.




Figure 3. EIC of select compounds at 0.01 ppm.

Initial Calibration

Calibration solutions of 0.01, 0.02, 0.05, 0.1, 0.2, 0.5, 1, 2, 5 and 10 ppm were prepared for all the analytes. The %RSD for each analyte should be less than 30% in order to use the mean response factors or linear regression fitting for calculating results. Quadratic calculations may not be used in EPA 525.2.




Figure 4. Calibration of benzo(a)pyrene from 0.01 to 10 ppm.

All compounds meet the method QC criteria. They showed excellent calibration coefficient and relative standard deviation at concentrations ranging from 0.01 to 10 ppm (Table 1), even for late eluting compounds such as benzo[a] pyrene (Figure 4). The average r^2 and %RSD of all 25 regulated compounds are 0.9975 and 14.58%, respectively.

Table 1. Calibration data of 25 regulated compounds.

| Compound Name | Correlation Coefficient (r^2) | Average RF | %RSD |
|-----------------------------|-----------------------------------|------------|-------|
| Hexachlorocyclopentadiene | 0.9988 | 0.6588 | 11.06 |
| 2,6-Dinitrotoluene | 0.999 | 0.7087 | 12.97 |
| 2,4-Dinitrotoluene | 0.9997 | 0.6845 | 16.34 |
| Propachlor | 0.9953 | 2.1286 | 12.65 |
| Hexachlorobenzene | 1.0000 | 0.5952 | 6.71 |
| Simazine | 0.9925 | 0.3924 | 11.08 |
| Atrazine | 0.9957 | 0.6690 | 16.17 |
| Lindane | 0.9999 | 1.7157 | 10.46 |
| Metribuzin | 0.9979 | 0.5000 | 13.6 |
| Alachlor | 0.9996 | 0.8753 | 9.53 |
| Heptachlor | 0.9985 | 0.7339 | 8.98 |
| Metolachlor | 0.9995 | 1.2759 | 8.63 |
| Cyanazine | 0.9958 | 0.5948 | 13.48 |
| Aldrin-R | 0.9941 | 0.9424 | 14.94 |
| Heptachlor epoxide | 0.9980 | 0.3972 | 8.72 |
| Butachlor | 0.9990 | 1.1007 | 11.55 |
| gamma-Chlordane | 0.9982 | 0.5831 | 14.45 |
| alpha-Chlordane | 0.9997 | 0.8379 | 9.02 |
| trans Nonachlor | 0.9997 | 0.2656 | 9.15 |
| Dieldrin | 0.9995 | 1.3265 | 8.92 |
| Endrin | 0.9996 | 0.1338 | 46.59 |
| bis(2-Ethylhexyl) adipate | 0.9963 | 3.2848 | 17.38 |
| Methoxychlor | 0.9924 | 1.3232 | 29.65 |
| bis(2-Ethylhexyl) phthalate | 0.9912 | 1.7660 | 29.84 |
| Benzo[a]pyrene | 0.9972 | 1.9991 | 12.6 |
| Average | 0.9975 | 1.0197 | 14.58 |

Method detection limits (MDLs) of these 25 compounds were calculated based on the standard deviation of seven replicates at 0.01 ppm multiplied by Student's t at 99% confidence level (Table 2).

Table 2. MDL replicates at 0.01 ppm.

| Compound Name | Average Amount (ppm) | %RSD | MDL (ppm) |
|------------------------------|----------------------|-------|-----------|
| Hexachlorocyclopentadiene | 0.0100 | 7.52 | 0.002372 |
| 2,6-Dinitrotoluene | 0.0100 | 14.48 | 0.004568 |
| 2,4-Dinitrotoluene | 0.0088 | 14.78 | 0.004087 |
| Propachlor | 0.0145 | 3.94 | 0.00179 |
| Hexachlorobenzene | 0.0126 | 21.21 | 0.008421 |
| Simazine | 0.0170 | 7.48 | 0.003996 |
| Atrazine | 0.0166 | 6.47 | 0.003372 |
| Lindane | 0.0133 | 6.07 | 0.002537 |
| Metribuzin | 0.0123 | 12.54 | 0.004859 |
| Alachlor | 0.0087 | 5.19 | 0.001418 |
| Heptachlor | 0.0115 | 13.49 | 0.004862 |
| Metolachlor | 0.0100 | 8.03 | 0.002522 |
| Cyanazine | 0.0100 | 26.48 | 0.008337 |
| Aldrin-R | 0.0098 | 11.45 | 0.003514 |
| Heptachlor epoxide | 0.0087 | 10.85 | 0.002957 |
| Butachlor | 0.0110 | 15.95 | 0.005505 |
| gamma-Chlordane | 0.0076 | 8.74 | 0.002097 |
| alpha-Chlordane | 0.0082 | 6.76 | 0.001739 |
| trans Nonachlor | 0.0091 | 13.47 | 0.003864 |
| Dieldrin | 0.0113 | 16.10 | 0.005729 |
| Endrin | 0.0108 | 28.40 | 0.009107 |
| bis(2-Ethylhexyl) adipate | 0.0116 | 9.07 | 0.003296 |
| Methoxychlor | 0.0048 | 14.22 | 0.002135 |
| bis(2-Ethylhexyl) phthalate* | 0.0267 | 9.30 | 0.007796 |
| Benzo[a]pyrene | 0.0098 | 11.26 | 0.003484 |
| Average | 0.0100 | 12.13 | 0.0042 |

*bis(2-Ethylhexyl) phthalate is at 0.05 ppm

Conclusion

The Varian 240-MS showed excellent sensitivity and linearity for all 25 regulated analytes. The 240-MS has the best full scan sensitivity of all GC/MS products on the market. V:Results GC/MS software allowed for automatic calculation of the MDLs. The average MDL of all 25 analytes is 3.5 ppb, significantly lower than that required in the EPA method.

These data represent typical results.

For further information, contact your local Varian Sales Office.

FactorFour, V:Results, Varian and the Varian Logo are trademarks or registered trademarks of Varian, Inc. in the U.S. and other countries.

© 2010 Varian, Inc.

Application Note SI-02374

Varian, Inc.

www.varianinc.com

North America: 800.926.3000 – 925.939.2400

Europe: *The Netherlands*: 31.118.67.1000

Asia Pacific: *Australia*: 613.9560.7133

Latin America: *Brazil*: 55.11.3238.0400



VARIAN