



Styrene/Isoprene Block Copolymers with Agilent PLgel MIXED-D and GPC/SEC

Application Note

Materials Testing and Research, Polymers

Authors

Greg Saunders and Ben MacCreath
Agilent Technologies (UK) Ltd
Essex Rd
Church Stretton
SY6 6AX
UK

Introduction

Unlike random copolymers where two or more monomers are included in a polymer chain in a statistically random manner, block copolymers contain isolated blocks of what are essentially homopolymers linked at one or more terminal positions. This regulated structure produces materials with novel properties not observed in random copolymer analogues, typically associated with phase behaviour at interfaces. Such copolymers have been exploited for their surfactant characteristics and for their propensity to form isolated domains in thin films. The simplest block copolymers are AB diblocks consisting of two chains of different polymers (A and B) joined at one terminal. Triblock polymers contain three separate polymer chains connected in series, typically composed of two (ABA triblocks) or three (ABC) different types of polymer chain. Producing block copolymers from blocks with widely differing chemistries results in materials with very unusual properties.

A family of block copolymers of styrene (S) and isoprene (I) have been developed for their excellent adhesive properties. In these materials, SI diblock and SIS triblock copolymers are blended to give a product with the required viscosity and modulus for particular applications.

Gel permeation chromatography (GPC) can be used to estimate the relative proportion of the two components in an SI/SIS product.



Agilent Technologies

Analysis of Styrene/Isoprene Block Copolymers

In this example, the relative styrene chain lengths in the diblock and triblock materials are the same and so the diblock will have a reduced molecular weight and therefore size in solution compared to the triblock.

Using high efficiency GPC columns, the two components can be resolved based on this size difference. Figure 1 shows a chromatogram of a commercial SI/SIS blended product. The SI and SIS components have been resolved allowing quantification to take place. In this case, the SI content is approximately 42%, typical of a material that finds application in reducing the modulus of polyolefins.

| | |
|-----------|--|
| Sample | Styrene isoprene block copolymer blend |
| Columns | 2 × Agilent PLgel 5 µm MIXED-D, 300 × 7.5 mm (p/n PL1110-6504) |
| Eluent | THF (stabilized) |
| Flow rate | 1.0 mL/min |
| Inj Vol | 100 µL |
| Detector | RI |
| System | Agilent PL-GPC 50 |




Figure 1. Chromatogram of a styrene isoprene di- and triblock copolymer blend on an Agilent PLgel 5 µm MIXED-D, column.

Conclusions

High efficiency Agilent PLgel MIXED-D columns resolve components of an isoprene/styrene block copolymer based on their size in solution. This allows quantification to take place. Typically, an increase in the proportion of diblock in the blend reduces the viscosity and modulus of the product, a feature that is particularly suited to pressure-sensitive label applications, for example. In this case, the styrene/isoprene content is approximately 42%, typical of a material used to reduce the modulus of polyolefins.

For More Information

These data represent typical results. For more information on our products and services, visit our Web site at www.agilent.com.

www.agilent.com

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2011
Printed in the USA
May 18, 2011
5990-8228EN



Agilent Technologies