


# Analysis of Phenols in Water with Star 9080 EC Detector

# LC

Varian Application Note  
Number 19

Fred Klink  
Varian Chromatography Systems

**Key Words:** Star 9080, phenols, environmental




**Analysis of a standard mixture of priority pollutants as listed by the US Environmental Protection Agency (EPA).**

The mixture contains: 2,4-dinitrophenol (1), phenol (2), 4-nitrophenol (3), 2-methyl-4,6-dinitrophenol (4), 2-chlorophenol (5), 2-nitrophenol (6), 2,4-dimethylphenol (7), 4-chloro-3-methylphenol (8), 2,4-dichlorophenol (9), and 2,4,6-trichlorophenol (10). Inset: polluted river water containing 8 ppb phenol.

<b>Detector</b>	Varian Star 9080 Amperometric Electrochemical Detector
<b>Column</b>	ODS, 100 x 4.6 mm, 3 $\mu$
<b>Flow rate</b>	1.5 mL/min
<b>Mobile phase</b>	50 mM HAC/NaAc, pH 4.0, 35% acetonitrile
<b>Sample</b>	100 - 1000 nM phenols, 20 $\mu$ L injection
<b>Temperature</b>	30 $^{\circ}$ C
<b>Flowcell</b>	2.74 mm Glassy Carbon working electrode
<b>E-cell</b>	1200 mV (vs. Ag/AgCl)

NOTICE: This document contains references to Varian. Please note that Varian, Inc. is now part of Agilent Technologies. For more information, go to [www.agilent.com/chem](http://www.agilent.com/chem).

# On-column Concentration of Phenols



- A: 20  $\mu$ L injection of 200 ppb standard
- B: 2000  $\mu$ L injection of 2 ppb standard
- C: 2000  $\mu$ L injection of 0.5 ppb standard