



Analysis of Reacetylated Polyvinyl Alcohol

Application Note

Materials Testing & Research, Polymers

Authors

Greg Saunders, Ben MacCreath
Agilent Technologies, Inc.

Introduction

Polyvinyl alcohols (PVOH) are water soluble polymers, which can be analyzed by aqueous size exclusion chromatography (SEC). However, if reacetylated, the resultant polyvinyl acetate (PVAc) is rendered THF soluble, and so Agilent PLgel MIXED-B columns are the columns of choice.



Agilent Technologies

PLgel 10 μm MIXED-B columns are designed for high MW polymer analysis and demanding eluent conditions. The PLgel 10 μm MIXED-B spans a wide range of molecular weights, up to 10 million, with a linear calibration curve. It is particularly useful for molecular weight distributions where slightly higher than average MWs are encountered. The 10 μm particle size provides good resolution with relatively low pressures for enhanced lifetimes in demanding conditions.

Conditions

Columns: 2 x PLgel 10 μm MIXED-B,
300 x 7.5 mm (p/n PL1110-6100)

Eluent: THF

Flow Rate: 1.0 mL/min

Loading: 0.25%, 100 μL

Detection: RI




Figure 1. Analysis of reacetylated poly(vinyl) alcohol using PLgel MIXED-B columns

www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2010

Published in UK, October 12, 2010

SI-01779



Agilent Technologies