

CHIRANA CHEESE LIFT



Patient chairs CHIRANA CHEESE LIFT with SK 1-08 CHIRANA CHEESE LIFT with SK 1-01

Dentist tables

Dentist tables 5N on pantographic arm

Cart dentist tables 5N





System of instruments hoses leading

Upper delivery system - Whips



Hanging system




Control panels


Button-control panel



Touch screen control panel







- Spray and power set up
- Light for instrument on/off
- Chip blower function on/off
- Reverse of MM on/off
- Giromatic function on/off
- Mode control (by foot control or manual from keyboard)
- Disinfection programm activation
- Glass filling
- Spittoon rinsing
- Stop function
- OP light switch on/off
- Chair control

Additional tray for dentist table





Tray table under dentist table

1x silicone tray


1x stainless tray

2x silicone tray

2x stainless tray



Spittoon block



Hygiene:

Clean water bottle system
Disisnfection system
Both systems
DCI bottles or Chirana bottles





Assistant tables


EASY 2N – two positions



EASY 4N – four positions



Control panel on assistant table



- Glass filling function
- Spittoon rinsing
- Stop function
- OP light switch on/off
- Chair control
- Switch buttons for public water or bottle

Arms for assistan table


Arm with adjustable height




Simple short arm



Extended arm



Optional accessories for column




Arms of trays on column

Short arm 355 mm



Normal arm 725 mm




Arms of dentist table



Pantographic arm

1 x slicone or 2x silicone tray
1x stainless or 2x stainless tray
with or w/o handle
with or w/o cup holder
with or w/o flacons holders





OP lights

FARO ALYA LED

adjustable light intensity 3,000 - 50,000 lux



FARO MAIA LED

adjustable light intensity 3,000 - 35,000 lux



LED - A

adjustable light intensity 7,000 - 38,000 lux



Foot controls

Multifunctional foot control FCR1-7- WL



Multifunctional foot control FCR1-7



Multifunctional foot pedal



Foot control with spray



Foot switch (only istruments activation)



Patient chair SK1-01 - upholsteries

