

NOTICE: This document contains references to Varian. Please note that Varian, Inc. is now part of Agilent Technologies. For more information, go to [www.agilent.com/chem](http://www.agilent.com/chem).



## APPLICATION NOTE

Application 1935 - HPLC

### Flavonoids

#### Ginkgo Biloba

Column: Polaris™ 5 $\mu$  C18-A, 150 x 4.6mm  
Part nr.: A2000150X046  
Eluent: MeCN:H<sub>2</sub>O with 0.5% H<sub>3</sub>PO<sub>4</sub>, 50:50  
Flow: 1.0 ml/min  
Temp: 35°C  
Detection: UV @ 254nm  
Sample: 1. Isorhamnetin  
2. kaempferol

