

NOTICE: This document contains references to Varian.
Please note that Varian, Inc. is now part of Agilent
Technologies. For more information, go to
www.agilent.com/chem.



APPLICATION NOTE

Cysteine underivatized

Application 1997 LC

Column : Polaris 5 μ C18-A, 250 x 4.6mm
Part No.: A2000250X046
Eluent : A: 50 mM NH₄H₂PO₄ + 1 mM sodium heptane sulfonate, pH 2.5
B: MeCN A : B 95 : 5
Flow : 1.0 ml/min Detection: UV @ 210 nm

Peak identification

1. L-cysteine

