

NOTICE: This document contains references to Varian. Please note that Varian, Inc. is now part of Agilent Technologies. For more information, go to www.agilent.com/chem.



APPLICATION NOTE

Application 1930 - HPLC

Antibiotics

Potential LC/MS Method for Separation of Antibiotics

Column: Polaris™ 5μ C18-A, 50 x 2.0mm
Part nr.: A2000050X020
Eluent: A: 0.1% formic acid
B: MeCN+0.1% formic acid
Gradient: 0 to 50%B in 10 min
Flow: 0.8 ml/min
Detection: UV @ 215nm
Sample: 1. Cephaloridine
2. Penicillin G
3. Penicillin V

