



Tricyclic Antidepressants

Tricyclic Antidepressants

Column: Polaris C18A, 150 x 2mm, 3µm

Part nr.: A2001150X020

Eluent: A: ACN, B: 0,05% CH₃COOH

Gradient: 0min: 10 + 90%B

25 min: 20 + 80%B

30 min: 10 + 90%B

Flow: 0,3 ml/min

Detection: UV 270nm

Sample: 1. zeatin

2. zeatin riboside

3. kinetin

4. 2ip =dimethylaminopurin

5. BAP= 6benzyl-amino-purin

6. 2ipR= dimethylaminopurin riboside

7. BAPR=6benzyl-amino-purin riboside

