

NOTICE: This document contains references to Varian.  
Please note that Varian, Inc. is now part of Agilent  
Technologies. For more information, go to  
[www.agilent.com/chem](http://www.agilent.com/chem).



## APPLICATION NOTE

### Tryptic Digest of b-Lactoglobulin A

Application 1998 LC

Column : Polaris 5 $\mu$  C18-A, 250 x 4.6mm  
Part No.: A2000250X046  
Eluent : A: H<sub>2</sub>O + 0.1% TFA B : MeCN + 0.085% TFA  
Gradient : 5% to 45%B in 40 mins  
Flow : 1.0 ml/min  
Detection : UV @ 215 nm

Sample : Tryptic digest of bovine b-lactoglobulin A

