

NOTICE: This document contains references to Varian.
Please note that Varian, Inc. is now part of Agilent
Technologies. For more information, go to
www.agilent.com/chem.



APPLICATION NOTE

Vitamin B6 complex

Application 2236 LC

Separation of Vitamin B6 complex on Pursuit Diphenyl

Column : Pursuit Diphenyl
Dimensions : 250 x 4.6 mm, 5 μ , Part No. A3040250x046
Mobile Phase : 5 mM Ammonium formate, pH 2.65:CH₃CN - 95:5
Flow Rate : 0.5 mL/min
Temperature : 40^oC
Detection : 254 nm

1. Pyridoxamine
2. Pyridoxal
3. Pyridoxine (or Pyridoxol or Vitamin B6)

